

PROSTIVA[®] RF THERAPY



Price List

Effective October 1, 2007

Ordering Information

To Place an Order

Call Medtronic Neuromodulation Customer Service.

Toll-Free: 1-888-638-7627 Fax: 1-800-930-9361

Customer Service representatives are on duty to take your order between 8 a.m. and 5 p.m., Central Standard Time. For faster service when calling, please have your account number available.

Payment Terms

Net 30 days from date of invoice. Prices are subject to change without notice.

Shipping

PROSTIVA® RF Systems are shipped Freight on Board (F.O.B.) shipping point, where current pricing includes Federal standard air or equivalent. Please allow five (5) days for delivery. Rush (air freight) orders are shipped F.O.B. shipping point. Vendor absorbs all shipping costs for returns generated by vendor shipping error, including freight charges both ways.

Credit and Returns

Full credit will be given for unopened, undamaged, and unmarked packages returned within 90 days from date of invoice. RF Generators that have been opened or used may be returned within 90 days with a 25% restocking fee. Please call Medtronic Neuromodulation Customer Service for a **RGR#** prior to returning the product and be sure the RGR# is visible on the outside of the package. Returns due to overstock, restock, or undeliverable may be sent to:

Medtronic, Inc.
Attention: Returns
2240 Woodale Drive
Mounds View, MN 55112

Warning

The products described in this price list are intended for use only with Medtronic components. The use of other manufacturers' components with these products may result in damage to Medtronic components, less than adequate performance, or increased risks to the patient.

PROSTIVA RF Therapy for Benign Prostatic Hyperplasia (BPH)

PROSTIVA RF Therapy Product/Component Overview.....	2
PROSTIVA RF System Components.....	3
Index	

The following are registered trademarks of Medtronic, Inc., and its affiliated companies:

Medtronic®, PROSTIVA® RF

PROSTIVA RF Therapy

PLEASE READ BEFORE ORDERING

When ordering, please consider all the components needed for a complete system for each patient procedure. Refer to specific model numbers in this catalog when placing your order.

PROSTIVA RF Therapy Overview

PROSTIVA RF System

A PROSTIVA RF system consists of each of the items that are contained in the 8931KIT0/15 PROSTIVA RF Complete Kit [PROSTIVA RF Generator, Rod Lens Telescope, Power Cord, PROSTIVA RF System User Guide and a standard one-year limited warranty agreement], plus a PROSTIVA Disposable Handpiece Kit.



8930001 PROSTIVA RF Generator



8929001 PROSTIVA Integrated Disposable Handpiece Kit

PROSTIVA RF Therapy System Components

ORDERING NUMBER	DESCRIPTION	
893KIT0/15	<p>PROSTIVA RF Generator Kit with 0° or 15° Scope</p> <p>Kit includes the following items: (1) Model 8930001 PROSTIVA RF Therapy Generator, (1) Power Cord, (1) PROSTIVA RF System User Guide, (1) Model 8099TU 0°/15° Rod Lens Telescope, and (1) standard one-year limited warranty.</p>	
8930001	<p>PROSTIVA RF Generator</p> <p>Includes the following items: (1) Model 8930001 PROSTIVA RF Therapy Generator, (1) Power Cord, (1) PROSTIVA RF System User Guide, and (1) standard one-year limited warranty.</p>	
8099TU/15	<p>Rod Lens Telescope</p> <p>Compatible with all medical-grade video camera systems and light sources, and can be used with direct vision. The telescope is a rigid optic telescope with a 0°/15° viewing angle. Designed for reuse. Autoclave Sterilization. Content: (1) Eyepiece with Glass Ocular Lens, (1) Fiber Optic Light Cable Connection with screw-on type cable connector to accept Circon ACMI - Storz - Wolf - Olympus light cables, (1) Stainless Steel Shaft containing rod lens system and objective lens, literature.</p>	
8929001	<p>PROSTIVA Integrated Disposable Handpiece Kit</p> <p>Single-use disposable handpiece designed for use only with the Model 8930001 PROSTIVA RF Generator System. Sterile contents include: (1) Single-use Disposable Handpiece designed with integrated RF cable and RF activation switch, (1) Return Electrode, and (1) Irrigation System. Non-sterile Content: Literature.</p>	<p>1-5 6-12 13+</p>
60883	<p>Foot Switch (Optional Accessory)</p> <p>Compatible only with the 8930001 PROSTIVA RF Generator. Includes: (1) Model 60883 Foot Switch, (1) User Guide, and (1) standard one-year limited warranty.</p>	
8919 COMING SOON	<p>PROSTIVA System Two-Year Extended Limited Warranty</p> <p>PROSTIVA System two-year extended limited warranty covers (1) RF Generator, (1) Rod Lens Telescope, and (1) Optional PROSTIVA Foot Switch (if purchased).</p>	

Index

Ordering Number	Description	Price	Page
8099TU/15	Rod Lens Telescope		3
8919	PROSTIVA System Two-Year Extended Limited Warranty		3
8931KIT0/15	PROSTIVA RF Generator Kit with 0° or 15° Scope		3
60883	Foot Switch		3
8929001	PROSTIVA Integrated Disposable Handpiece Kit		3
8930001	PROSTIVA RF Generator		3

United States of America

Medtronic Neuromodulation
710 Medtronic Parkway
Minneapolis, MN 55432-5604
USA
Internet: www.medtronic.com
Tel. 763-505-5000
Fax 763-505-1000
Toll-free 1-800-328-0810

Europe

Medtronic International Trading Sàrl
Route du Molliau 31
Case Postale
CH-1131 Tolochenaz
Switzerland
Tel. +41-21-802-7000
Fax +41-21-802-7900

Asia-Pacific

Medtronic International, Ltd.
Suite 1602 16/F
Manulife Plaza
The Lee Gardens, 33 Hysan Avenue
Causeway Bay
Hong Kong
Tel. 852-2891-4456
Fax 852-2891-6830

Australia

Medtronic Australasia Pty. Ltd.
97 Waterloo Road
North Ryde
NSW 2113
Australia
Tel. +61-2-9857-9000
Fax +61-2-9878-5100
www.medtronicneuro.com.au

Canada

Medtronic of Canada Ltd.
6733 Kitimat Road
Mississauga, Ontario L5N 1W3
Canada
Tel. 1-905-826-6020
Fax 1-905-826-6620

UC200300919b EN N15393b
© Medtronic, Inc. 2007
All Rights Reserved
Printed in USA

**Medtronic***Alleviating Pain · Restoring Health · Extending Life*